

ContentTag	Original Content	Latest Content	Comments	Date Modified	-	KEY
ST1	End User License Agreement		Screen Title	07/17/18 1:49 PM		T# = Text
ST2	Put Your Implant Into Discoverable Mode		Screen Title	07/17/18 1:49 PM		B# = Button
ST3	Discover & Pair		Screen Title	07/17/18 1:49 PM		L# = Label
ST4	No HF10 Implants Discovered	HF10 Implant not discovered	Screen Title	07/31/18 11:10 AM		I# = Illustration
ST5	Unable to Pair		Screen Title	07/17/18 1:49 PM		A# = Animation
ST6	Name	What is your name?	Screen Title	07/31/18 11:39 AM		N# = Navigation
ST7	Confirm Pain Areas	Confirm Your Pain Areas	Screen Title	07/31/18 11:39 AM		F# = Form Field
ST8	I Want To Be Free To...		Screen Title	07/17/18 1:49 PM		ST# = Screen Title
ST9	Notification Settings		Screen Title	07/19/18 7:47 PM		In Progress
ST10	App Security		Screen Title	07/19/18 7:47 PM		
ST11	Change HF10 Settings		Screen Title	08/02/18 1:42 PM		
ST12	My Progress		No longer needed with new design.	09/11/18 11:00 PM		
ST13	Favorites		Section Title	07/18/18 2:23 PM		
ST14	Menu		Section Title	07/18/18 2:23 PM		
ST15	HF10 Optimization Survey		Section Title	08/02/18 9:21 AM		
ST16	HF10 Pain Survey		Section Title	08/02/18 9:21 AM		
ST17	Select Program		Workflow Title	09/03/18 6:35 PM		
B1	Accept		Button Label to agree with EULA	07/17/18 12:51 PM		
B2	Decline		Button label to disagree with EULA	07/17/18 12:51 PM		
B3	Return to Agreement			07/17/18 12:51 PM		
B4	Re-check			07/17/18 12:51 PM		
B5	?		Button Label for Support/Help	07/31/18 12:35 PM		
B6	Pair		Button Label	07/17/18 12:51 PM		
B7	Start Over		Button Label	07/17/18 12:51 PM		
B8	Cancel		Button Label	07/17/18 12:51 PM		
B9	Continue		Button Label	07/17/18 12:51 PM		
B10	Next		Button Label	07/17/18 12:51 PM		
B11	Setup Later in Settings	Set up Later in Settings	Button Label	07/31/18 11:57 AM		
B12	(X within Circle)		Asset can be found in iOS/Android UI Kit in Zeplin	07/19/18 7:55 PM		
B13	Done		Button Label when in Edit Field state with OSK visible.	07/17/18 12:51 PM		
B14	More About Pain Areas		Contextual Help Button Label	07/17/18 12:51 PM		
B15	dance		Button Label	07/17/18 12:51 PM		
B16	take a trip		Button Label	07/17/18 12:51 PM		
B17	ride a bike		Button Label	07/17/18 12:51 PM		
B18	play catch		Button Label	07/17/18 12:51 PM		
B19	go out on the town		Button Label	07/17/18 12:51 PM		
B20	+ add yours		Begins B20 Custom Goal Sub-workflow	07/17/18 12:51 PM		
B21	Save		Button Label for individual field being edited.	07/17/18 12:51 PM		
B22	Allow Push Notifications			07/17/18 12:51 PM		
B23	Don't Allow			07/17/18 12:51 PM		
B24	OK			07/17/18 12:51 PM		
B25	Setup PIN		Button Label	07/17/18 12:51 PM		
B26	Enable Touch ID		Button Label	07/17/18 12:51 PM		
B27	HF10 Logo inside a circle		Tapping expands/collapses Control Panel	07/18/18 11:05 AM		
B28	Change HF10 Settings			07/18/18 11:05 AM		
B29	Done			07/18/18 12:23 PM		
B30	Dismiss			07/19/18 3:34 PM		
B31	Exit			07/27/18 4:17 PM		
B32	Assessing Pain Relief			07/27/18 4:27 PM		
B33	Apply Recommendation			07/30/18 11:55 AM		
B34	Stay on Current Settings			07/30/18 11:55 AM		
B35	My Progress			07/30/18 11:55 AM		
B36	Home			07/30/18 11:55 AM		
B37	Tips & Guidance			07/30/18 12:06 PM		
B38	Return			07/30/18 12:34 PM		
B39	Exit Survey			07/30/18 12:34 PM		
B40	Assessing Pain Scores			07/30/18 1:16 PM		
B41	Activity with HF10			07/30/18 1:16 PM		
B42	Pain Medication with HF10			07/30/18 1:16 PM		
B43	Take a Survey			10/10/18 10:58 PM		
B44	Back			07/31/18 10:42 AM		
B45	Body			08/28/18 10:38 AM		
B46	Graph			08/28/18 10:38 AM		
B47	Program History			08/28/18 10:38 AM		
B48	Activate			08/28/18 4:04 PM		
B49	Favorite			08/28/18 4:04 PM		
B50	More About My Progress			08/28/18 4:07 PM		
B51	Delete			09/03/18 5:17 PM		
B52	Edit			09/03/18 5:30 PM		
B53	Add			09/03/18 6:35 PM		
B54	Deactivate Favorite			10/10/18 10:58 PM		
	Begin Troubleshooting			10/10/18 10:58 PM		
L1	Acceptance of Terms		Alert Title for user who declines EULA - Nevro needs to review and determine if it is necessary	07/17/18 12:51 PM		
L2	Phone Incompatible			07/17/18 12:51 PM		
L3	HF10 Implant #X [where X is serial number]	HF10 Implant #X [where X is serial number]		07/31/18 11:16 AM		
L4	Why Allow Notifications?		Help Topic Title	07/17/18 12:51 PM		
L5	Enable Touch ID or PIN		Label Text	07/17/18 12:51 PM		
L6	Enter a 4-digit PIN		Instructional Text	07/17/18 12:51 PM		
L7	Verify your new PIN		Instructional Text	07/17/18 12:51 PM		
L8	PIN Saved		Confirmation Label	07/17/18 12:51 PM		
L9	[Warning Icon] PINs did not match. Try again.		Error text with instructions	07/17/18 12:51 PM		
L10	Congratulations!		Confirmation Label	07/17/18 12:51 PM		
L11	Last Data Sync		Removed from old control panel and re-purposing on Menu Screen: App Info	10/08/18 2:05 PM		
L12	Program			07/18/18 11:20 AM		
L13	Level			07/18/18 11:20 AM		
L14	Implant Battery		Status Indicator Label; because this is no longer visually rendered in the UI, still keeping it since it should be spec'd as label for screen readers (accessibility).	10/08/18 3:07 PM		

L15	BLE Connection	Bluetooth Connection	Status Indicator Label; because this is no longer visually rendered in the UI, still keeping it since it should be spec'd as label for screen readers (accessibility).	10/08/18 3:07 PM		
L16	Default Level		Label for disabled Level Adjustment	07/18/18 12:15 PM		
L17	HF10 is Currently Off		Alert Title for user changing program while HF10 is off.	07/18/18 12:37 PM		
L18	HF10 Battery Low			07/19/18 3:34 PM		
L19	HF10 Battery Critically Low		GH: Removed 10/5 since 25, 10 and 5% warnings are being modified to be "Low" and "Depleted"	10/05/18 12:53 PM		
L20	HF10 Battery Depleted			07/19/18 3:34 PM		
L21	HF10 Battery Fully Charged			07/19/18 3:34 PM		
L22	HF10 BLE Connection Problem			07/19/18 3:46 PM		
L23	Overall Pain Relief			07/27/18 4:17 PM		
L24	No Relief			07/31/18 11:05 AM		
L25	Fair			07/27/18 4:53 PM		
L26	Good			07/31/18 11:05 AM		
L27	Very Good			07/31/18 11:05 AM		
L28	Excellent			07/31/18 11:05 AM		
L29				07/31/18 11:05 AM		
L30				07/31/18 11:05 AM		
L31				07/31/18 11:05 AM		
L32				07/31/18 11:05 AM		
L33	Pain Scores (n of x)		Where n = 1, 2, 3, 4 or 5 out of x total number of pain areas that are being targeted by HF10. Minimum of 1, maximum of 5)	07/30/18 2:07 AM		
L34	Activity Level			07/30/18 11:08 AM		
L35	Decreased			07/30/18 11:12 AM		
L36	Stayed the Same			07/30/18 11:12 AM		
L37	Increased			07/30/18 11:12 AM		
L38	Pain Medication			07/30/18 11:20 AM		
L39	HF10 Recommendation			07/30/18 11:53 AM		
L40	Thank You!			07/30/18 12:06 PM		
L41	Exit Survey?			07/30/18 12:34 PM		
L42	Testing [X]		Where [X] = (Program Number + Level) for PTA or (Setting Name) for Pathfinder. e.g. P2 L1 P5 L7 JD15 (where JD are the patient's initials under Pathfinder).	10/10/18 10:48 PM		
L43	Survey Reminder	Reminder		08/02/18 8:29 AM		
L44	Time to Optimize Your HF10			07/30/18 6:04 PM		
L45	Time to Log Your Progress			07/30/18 6:04 PM		
L46	Current Relief			07/30/18 6:04 PM		
L47	Start Date			07/30/18 6:04 PM		
L48	Time			07/30/18 6:04 PM		
L49	Settings			07/30/18 6:04 PM		
L50	Pain Scores	Pain Score per Area		10/05/18 6:33 PM		
L51	Overall Relief Score			08/28/18 10:53 AM		
L52	No Data			08/28/18 3:53 PM		
L53	Pain Area			08/28/18 3:57 PM		
L54	Score			08/28/18 3:57 PM		
L55	Baseline			08/28/18 3:57 PM		
L56	HF10 Off			08/28/18 4:37 PM		
L57	HF10 On			08/28/18 4:37 PM		
L58	Active Favorite			09/03/18 5:40 PM		
L59	Add Favorite			09/03/18 6:03 PM		
L60	Add Current Program			09/03/18 6:03 PM		
L61	Add from Program History			09/03/18 6:03 PM		
L62	Edit Favorite			09/03/18 6:44 PM		
L63	Delete Favorite			09/03/18 7:01 PM		
L64	Unable to Add Favorite			09/03/18 7:01 PM		
L65	Setting			10/08/18 1:40 PM		
L66	"[X]" Favorite		Where [X] = Favorite Name that's currently running for either Pathfinder or PTA is displayed in quotes, followed by the word Favorite. e.g., "Walking" Favorite "Running" Favorite "Sleeping" Favorite	10/09/18 1:37 PM		
L67	HF10 is currently ON.			10/09/18 12:46 PM		
L68	Favorite in Use			10/10/18 10:58 PM		
L69	[X] Optimized		Where [X] = Pathfinder Program Name. e.g. "JD2 Optimized"	10/10/18 10:58 PM		
L70	Connection Lost			10/10/18 10:58 PM		
L71	Status		Label at top of home screen cards	10/10/18 10:58 PM		
T1	EULA Agreement Copy		Currently with Nevro Legal	07/17/18 12:51 PM		
T2	To use this app for controlling your HF10 Implant, you must agree to the End User License Agreement on the screen. If you are certain you wish to decline acceptance of these terms, please quit and uninstall the app from your phone.		Text for user who declines EULA - Nevro needs to review and determine if it is necessary	07/17/18 12:51 PM		
T3	Your phone isn't compatible with HF10. You can re-check your phone or seek Support / Help.			07/17/18 12:51 PM		
T4	Scanning for HF10 Implant...			07/17/18 12:51 PM		
T5	Found HF10 Implant	HF10 Implant Found	Header	07/31/18 11:16 AM		
T6	Tap on your HF10 Implant below to pair, ensuring the serial number matches the one on the card given to you after surgery.	Confirm the serial number matches your patient ID, then tap the PAIR button below	Instructional Text	07/31/18 11:16 AM		
T7	Don't see your device? Tap below to start over.	Don't see your device?	Instructional Text	07/31/18 11:16 AM		
T8	Attempting to Pair...		System Status	07/17/18 12:51 PM		

T9	Congratulations!		Header	07/17/18 12:51 PM		
T10	Your HF10 Implant has been paired successfully.		System Status	07/17/18 12:51 PM		
T11	"Be sure your HF10 Implant is in "Discoverable" Mode. It will stay in this mode for two minutes. Tap below to start over and try again.		Instructional Text	07/17/18 12:51 PM		
T12	We were unable to pair with your HF10 Implant because it is no longer in "Discoverable" mode. It will stay in this mode for two minutes. Tap below to start over and try again.	We were unable to pair with your HF10 Implant because it is no longer in "Discoverable" mode. It will stay in this mode for two minutes.	Instructional Text	07/31/18 11:29 AM		
T13	Please provide your first name		Instructional Text	07/17/18 12:51 PM		
T14	mm/dd/yyyy, hh:mm AM/PM		Date and Time for Last Data Sync (L11)	10/08/18 2:05 PM		
T15	Your changes will be applied, but stimulation will remain OFF until you turn it ON.		Text for user making program level change while stim is off. GH: No longer relevant as of 10/5/2018.	10/05/18 12:30 PM		
T16	Slide to Turn Off		Text instructions for turning stimulation off. No longer relevant in new UI.	10/08/18 3:07 PM		
T17	Slide to Turn On		Text instructions for turning stimulation on. No longer relevant in new UI.	10/08/18 3:07 PM		
T18	Your HF10 is setup to target the pain area(s) above	Your HF10 is set up to target the pain area(s) above	Instructional Text	07/31/18 11:45 AM		
T19	Select activities you'd like to do.	Select activities you'd enjoy	Instructional Text	07/31/18 11:52 AM		
T20	Allow push notifications so your HF10 can notify you with low battery alerts and survey reminders.		Value prop statement for allowing notifications.	07/17/18 12:51 PM		
T21	To turn Notifications on later, go to Menu > Notification Settings and select the options you'd like. You can also enable them from your iOS Settings > Notifications and scrolling to HF10.		Help Topic Content	07/17/18 12:51 PM		
T22	To prevent accidental changes to your therapy settings, you can choose to setup Touch ID or PIN		Instruction Text	07/17/18 12:51 PM		
T23	Your PIN has been saved		Confirmation Text	07/17/18 12:51 PM		
T24	Your phone is now setup for use with your HF10 Implant		Confirmation Text	07/17/18 12:51 PM		
T25	Please recharge soon.			10/05/18 12:53 PM		
T26	Your HF10 battery is at 10%. Please recharge soon.		GH: Removed 10/5 since 25, 10 and 5% warnings are being modified to be "Low" and "Depleted"	10/08/18 3:11 PM		
T27	Please recharge now.			10/05/18 1:34 PM		
T28	Your HF10 battery is now fully charged.			07/19/18 3:46 PM		
T29	There is a problem connecting to your HF10.			07/19/18 3:43 PM		
T30	What is your overall pain relief today?			07/27/18 4:18 PM		
T31	What is your X pain score today?		The question is asked n number of times, with each time X equal to the respective pain area name from T32 array below.	07/30/18 2:23 AM		
T32[0—25]	left foot, right foot, bilateral feet, left leg, right leg, bilateral leg, perineum, left buttock, right buttock, bilateral buttock, lower back, mid back, upper back, neck, left hand, right hand, left arm, right arm, bilateral arm, left shoulder, right shoulder, bilateral shoulder, groin, abdomen, chest, face		An array containing all 26 distinct pain areas names for HF10.	07/30/18 2:23 AM		
T33	With HF10, your activity level has:			07/30/18 11:12 AM		
T34	With HF10, your pain medication has:			07/30/18 11:20 AM		
T35	Based on your answers, there is a recommended setting you should try.			07/30/18 11:53 AM		
T36	Based on your answers, it's recommended you stay on our current settings. Your responses are correlated with data from your HF10 Implant, and you can view this in My Progress.			07/30/18 11:53 AM		
T37	Based on your answers, it's recommended you test this program for another 48 hours. Your responses are correlated with data from your HF10 Implant, and you can view this in My Progress.			07/30/18 11:53 AM		
T38	Your responses are correlated with data from your HF10 Implant, and you can view this in My Progress.			07/30/18 12:06 PM		
T39	To maximize the benefit you get from HF10, we encourage you to check out more information in Tips & Guidance.			07/30/18 12:06 PM		
T40	Are you sure you wish to exit the survey? Your responses won't be saved.			07/30/18 12:34 PM		
T41	48 hours is recommended when trying out new program settings.	Staying with each program for 48 hours is recommended when trying out new program settings.		07/31/18 12:17 PM		
T42	Take the short survey so we can provide recommendations to optimize your HF10.			07/30/18 6:15 PM		
T43	Take the short survey so you can log your progress with HF10.			07/30/18 6:15 PM		
T44	We encourage you to continue testing the current program for another X hours.			08/28/18 11:34 AM		
T45	Turn on HF10 and begin taking surveys to see this data mapped above.			08/28/18 3:53 PM		
T46	Turn on HF10 and begin taking surveys to see your results graphed over time.			08/28/18 4:37 PM		
T47	No survey was taken, so we don't know what % pain relief you experienced during this program.			08/28/18 4:37 PM		
T48	You can save up to 10 favorites			09/03/18 5:12 PM		

T49	You can save [X] more favorites		Where X is the number of remaining favorite slots (10 total) that are available.	09/03/18 5:14 PM		
T50	You can only save 10 favorites. Please delete one before adding another.			09/03/18 5:14 PM		
T51	Choose a custom name (12 character max) and pick a color to label your favorite.			09/03/18 6:08 PM		
T52	Are you sure you want to delete [X]?		Where X is the name of the Favorite the patient is wanting to delete.	09/03/18 7:01 PM		
T53	1. Unplug the charger		Frame 1 in Pairing Animation	10/07/18 6:41 PM		
T54	2. Press the power button until you hear a beeping sound		Frames 2 & 3 in Pairing Animation	10/07/18 6:41 PM		
T55	3. Place charging coil over implant site on body		Frame 4 in Pairing Animation	10/07/18 6:41 PM		
T56	4. When aligned, chirping will stop and battery indicator will light up		Frames 5 & 6 in Pairing Animation	10/07/18 6:41 PM		
T57	Tap to turn HF10 ON			10/09/18 12:41 PM		
T58	Your HF10 has been turned OFF.			10/09/18 12:46 PM		
T59	Your HF10 has been turned ON.			10/09/18 12:46 PM		
A1	Animation: How to Discover & Pair Your HF10 Implant		Looping animation depicting the steps for putting Implant into Discoverable Mode. See animation_discover-pair-[1-10] in Zeplin for animation specs	07/19/18 7:42 PM		
A2	Animation: HF10 Pairing Attempt		Animation Specs in Progress	07/19/18 7:42 PM		
I1	Illustration: HF10 with checkmark		Placeholder image, use illustration in Zeplin comp as working placeholder	07/19/18 8:04 PM		
I2	Illustration: HF10 with question mark icon		Placeholder image, use illustration in Zeplin comp as working placeholder	07/19/18 8:04 PM		
I3	Illustration: HF10 with warning icon		Placeholder image, use illustration in Zeplin comp as working placeholder	07/19/18 8:04 PM		
I4	Illustration: Generic type of activity depicted.		Placeholder image, use illustration in Zeplin comp as working placeholder	07/19/18 8:04 PM		
I5	Illustration: iOS Phone with HF10 Notification Banner. Text is: It's Time to Optimize your HF10. Take the Survey.		Placeholder image, use illustration in Zeplin comp as working placeholder. Illustration for what a notification looks like with HF10.	07/19/18 8:04 PM		
I6	Illustration: Why Allow Notifications		Placeholder image, use illustration in Zeplin comp as working placeholder	07/19/18 8:04 PM		
I7	Illustration: Touch ID or other graphic representing authentication on iOS / Android.		Placeholder image, use illustration in Zeplin comp as working placeholder	07/19/18 8:04 PM		
I8	Illustration: Trophy cup		Placeholder image, use illustration in Zeplin comp as working placeholder	07/19/18 8:04 PM		
I9	Illustration: Posterior view of gender-neutral silhouette of human body upon which pain area overlays will be placed.		See iOS/Android UI Kit in Zeplin	07/19/18 7:52 PM		
I10	Illustration: Pain area overlay that will be used to visually indicate a pain area.		See iOS/Android UI Kit in Zeplin	07/19/18 7:52 PM		
I11	Illustration I9 + Anterior view of same silhouette and placed next to it when pain areas from both front AND back are being targeted.		See iOS/Android UI Kit in Zeplin	07/19/18 7:52 PM		
I12	Illustration: Activity Level Lifestyle Image for Survey Question		-	07/30/18 11:58 AM		
I13	Illustration: Pain Medication Lifestyle Image for Survey Question		-	07/30/18 11:58 AM		
I14	Illustration: Pathfinder Result 1		-	07/30/18 12:28 PM		
I15	Illustration: Pathfinder Result 2		-	07/30/18 12:28 PM		
I16	Swipe Left-Right Gesture		See iOS/Android UI Kit in Zeplin	07/19/18 7:55 PM		
I17	Program Adjustment		See iOS/Android UI Kit in Zeplin, Contains #s 1—5 and, representing each of the five programs programmed for the patient.	07/19/18 7:55 PM		
I18	Illustration: Pathfinder Result 3			07/30/18 12:28 PM		
I19	Illustration: PTA Result 3			07/30/18 12:28 PM		
I20	Illustration: PTA Result 3			07/30/18 12:34 PM		
I21	Illustration: Pathfinder Result 3			07/31/18 10:41 AM		
F1	Form Field		Empty by default	07/19/18 7:55 PM		
F2	Form Field		Appears after tapping B20	07/19/18 7:55 PM		
N1	Home		App Navigation	07/18/18 3:00 PM		
N2	My Progress		App Navigation	07/18/18 3:00 PM		
N3	Favorites		App Navigation	07/18/18 3:00 PM		
N4	Menu		App Navigation	07/18/18 3:00 PM		