

Content Tag	Original Content	Latest Content	Comments	Status	Date Modified	
ST1	End User License Agreement		Screen Title	Ready for review	07/17/18 1:49 PM	
ST2	Put Your Implant Into Discoverable Mode		Screen Title	Approved	07/17/18 1:49 PM	Color Key
ST3	Discover & Pair		Screen Title	Approved	07/17/18 1:49 PM	
ST4	No HF10 Implants Discovered	HF10 Implant not discovered	Screen Title	Ready for review	07/31/18 11:10 AM	ST# = Screen Title
ST5	Unable to Pair		Screen Title	Ready for review	07/17/18 1:49 PM	B# = Button
ST6	Name	What is your name?	Screen Title	Ready for review	07/31/18 11:39 AM	L# = Label
ST7	Confirm Pain Areas	Confirm Your Pain Areas	Screen Title	Ready for review	07/31/18 11:39 AM	T# = Text
ST8	I Want To Be Free To...		Screen Title	Ready for review	07/17/18 1:49 PM	A# = Animation
ST9	Notification Settings		Screen Title	Ready for review	07/19/18 7:47 PM	I# = Illustration
ST10	App Security		Screen Title	Ready for review	07/19/18 7:47 PM	F# = Form Field
ST11	Change HF10 Settings		Screen Title	Ready for review	08/02/18 1:42 PM	N# = Navigation
ST12	My Progress		No longer needed with new design.	Ready for review	09/11/18 11:00 PM	In Progress
ST13	Favorites		Section Title	Ready for review	07/18/18 2:23 PM	
ST14	Menu		Section Title	Ready for review	07/18/18 2:23 PM	
ST15	HF10 Optimization Survey		Section Title	Ready for review	08/02/18 9:21 AM	
ST16	HF10 Pain Survey		Section Title	Ready for review	08/02/18 9:21 AM	
ST17	Select Program		Workflow Title	Ready for review	09/03/18 6:35 PM	
B1	Accept		Button Label to agree with EULA	Drafting	07/17/18 12:51 PM	
B2	Decline		Button label to disagree with EULA	Drafting	07/17/18 12:51 PM	
B3	Return to Agreement			Drafting	07/17/18 12:51 PM	
B4	Re-check			Drafting	07/17/18 12:51 PM	
B5	?		Button Label for Support/Help	Drafting	07/31/18 12:35 PM	
B6	Pair		Button Label	Drafting	07/17/18 12:51 PM	
B7	Start Over		Button Label	Drafting	07/17/18 12:51 PM	
B8	Cancel		Button Label	Drafting	07/17/18 12:51 PM	
B9	Continue		Button Label	Drafting	07/17/18 12:51 PM	
B10	Next		Button Label	Drafting	07/17/18 12:51 PM	
B11	Setup Later in Settings	Set up Later in Settings	Button Label	Drafting	07/31/18 11:57 AM	
B12	(X within Circle)		Asset can be found in iOS/Android UI Kit in Zeplin	Drafting	07/19/18 7:55 PM	
B13	Done		Button Label when in Edit Field state with OSK visible.	Drafting	07/17/18 12:51 PM	
B14	More About Pain Areas		Contextual Help Button Label	Drafting	07/17/18 12:51 PM	
B15	dance		Button Label	Drafting	07/17/18 12:51 PM	
B16	take a trip		Button Label	Drafting	07/17/18 12:51 PM	
B17	ride a bike		Button Label	Drafting	07/17/18 12:51 PM	
B18	play catch		Button Label	Drafting	07/17/18 12:51 PM	
B19	go out on the town		Button Label	Drafting	07/17/18 12:51 PM	
B20	+ add yours		Begins B20 Custom Goal Subworkflow	Drafting	07/17/18 12:51 PM	
B21	Save		Button Label for individual field being edited.	Drafting	07/17/18 12:51 PM	
B22	Allow Push Notifications			Drafting	07/17/18 12:51 PM	
B23	Don't Allow			Drafting	07/17/18 12:51 PM	
B24	OK			Drafting	07/17/18 12:51 PM	
B25	Setup PIN		Button Label	Drafting	07/17/18 12:51 PM	
B26	Enable Touch ID		Button Label	Drafting	07/17/18 12:51 PM	
B27	HF10 Logo inside a circle		Tapping expands/collapses Control Panel	Drafting	07/18/18 11:05 AM	
B28	Change HF10 Settings			Drafting	07/18/18 11:05 AM	
B29	Done			Drafting	07/18/18 12:23 PM	
B29	Dismiss			Drafting	07/19/18 3:34 PM	
B30	Exit			Drafting	07/27/18 4:17 PM	
B31	Assessing Pain Relief			Drafting	07/27/18 4:27 PM	
B32	Apply Recommendation			Drafting	07/30/18 11:55 AM	
B33	Stay on Current Settings			Drafting	07/30/18 11:55 AM	
B34	My Progress			Drafting	07/30/18 11:55 AM	
B35	Home			Drafting	07/30/18 11:55 AM	
B36	Tips & Guidance			Drafting	07/30/18 12:06 PM	
B37	Return			Drafting	07/30/18 12:34 PM	
B38	Exit Survey			Drafting	07/30/18 12:34 PM	
B39	Assessing Pain Scores			Drafting	07/30/18 1:16 PM	
B40	Activity with HF10			Drafting	07/30/18 1:16 PM	
B41	Pain Medication with HF10			Drafting	07/30/18 1:16 PM	
B42	Take a Survey			Drafting	10/10/18 10:58 PM	
B43	Back			Drafting	07/31/18 10:42 AM	
B44	Body			Drafting	08/28/18 10:38 AM	
B45	Graph			Drafting	08/28/18 10:38 AM	
B46	Program History			Drafting	08/28/18 10:38 AM	
B47	Activate			Drafting	08/28/18 4:04 PM	
B48	Favorite			Drafting	08/28/18 4:04 PM	
B49	More About My Progress			Drafting	08/28/18 4:07 PM	
B50	Delete			Drafting	09/03/18 5:17 PM	
B51	Edit			Drafting	09/03/18 5:30 PM	

Content Tag	Original Content	Latest Content	Comments	Status	Date Modified
B52	Add			Drafting	09/03/18 6:35 PM
B53	Deactivate Favorite			Drafting	10/10/18 10:58 PM
B54	Begin Troubleshooting			Drafting	10/10/18 10:58 PM
L1	Acceptance of Terms		Alert Title for user who declines EULA - Nevro needs to review and determine if it is necessary	Drafting	07/17/18 12:51 PM
L2	Phone Incompatible			Drafting	07/17/18 12:51 PM
L3	HF10 Implant #X [where X is serial number]	HF10 Implant #X [where X is serial number]		Drafting	07/31/18 11:16 AM
L4	Why Allow Notifications?		Help Topic Title	Drafting	07/17/18 12:51 PM
L5	Enable Touch ID or PIN		Label Text	Drafting	07/17/18 12:51 PM
L6	Enter a 4-digit PIN		Instructional Text	Drafting	07/17/18 12:51 PM
L7	Verify your new PIN		Instructional Text	Drafting	07/17/18 12:51 PM
L8	PIN Saved		Confirmation Label	Drafting	07/17/18 12:51 PM
L9	[Warning Icon] PINs did not match. Try again.		Error text with instructions	Drafting	07/17/18 12:51 PM
L10	Congratulations!		Confirmation Label	Drafting	07/17/18 12:51 PM
L11	Last Data Sync		Removed from old control panel and re-purposing on Menu Screen: App Info	Drafting	10/08/18 2:05 PM
L12	Program			Drafting	07/18/18 11:20 AM
L13	Level			Drafting	07/18/18 11:20 AM
L14	Implant Battery Status		Indicator Label; because this is no longer visually rendered in the UI, still keeping it since it should be spec'd as label for screen readers (accessibility).	Drafting	10/08/18 3:07 PM
L15	BLE Connection	Bluetooth Connection	Status Indicator Label; because this is no longer visually rendered in the UI, still keeping it since it should be spec'd as label for screen readers (accessibility).	Drafting	10/08/18 3:07 PM
L16	Default Level		Label for disabled Level Adjustment	Drafting	07/18/18 12:15 PM
L17	HF10 is Currently Off		Alert Title for user changing program while HF10 is off.	Drafting	07/18/18 12:37 PM
L18	HF10 Battery Low			Drafting	07/19/18 3:34 PM
L19	HF10 Battery Critically Low		GH: Removed 10/5 since 25, 10 and 5% warnings are being modified to be "Low" and "Depleted"	Drafting	10/05/18 12:53 PM
L20	HF10 Battery Depleted			Drafting	07/19/18 3:34 PM
L21	HF10 Battery Fully Charged			Drafting	07/19/18 3:34 PM
L22	HF10 BLE Connection Problem			Drafting	07/19/18 3:46 PM
L23	Overall Pain Relief			Drafting	07/27/18 4:17 PM
L24	No Relief			Drafting	07/31/18 11:05 AM
L25	Fair			Drafting	07/27/18 4:53 PM
L26	Good			Drafting	07/31/18 11:05 AM
L27	Very Good			Drafting	07/31/18 11:05 AM
L28	Excellent			Drafting	07/31/18 11:05 AM
L29				Drafting	07/31/18 11:05 AM
L30				Drafting	07/31/18 11:05 AM
L31				Drafting	07/31/18 11:05 AM
L32				Drafting	07/31/18 11:05 AM
L33	Pain Scores (n of x)		Where n = 1, 2, 3, 4 or 5 out of x total number of pain areas that are being targeted by HF10. Minimum of 1, maximum of 5)	Drafting	07/30/18 2:07 AM
L34	Activity Level			Drafting	07/30/18 11:08 AM
L35	Decreased			Drafting	07/30/18 11:12 AM
L36	Stayed the Same			Drafting	07/30/18 11:12 AM
L37	Increased			Drafting	07/30/18 11:12 AM
L38	Pain Medication			Drafting	07/30/18 11:20 AM
L39	HF10 Recommendation			Drafting	07/30/18 11:53 AM
L40	Thank You!			Drafting	07/30/18 12:06 PM
L41	Exit Survey?			Drafting	07/30/18 12:34 PM
L42	Testing [X]		Where [X] = (Program Number + Level) for PTA or (Setting Name) for Pathfinder. e.g., P2 L1; P5 L7; JD15 (where JD are the patient's initials under Pathfinder).	Drafting	10/10/18 10:48 PM
L43	Survey Reminder		Reminder	Drafting	08/02/18 8:29 AM
L44	Time to Optimize Your HF10			Drafting	07/30/18 6:04 PM
L45	Time to Log Your Progress			Drafting	07/30/18 6:04 PM
L46	Current Relief			Drafting	07/30/18 6:04 PM
L47	Start Date			Drafting	07/30/18 6:04 PM
L48	Time			Drafting	07/30/18 6:04 PM
L49	Settings			Drafting	07/30/18 6:04 PM
L50	Pain Scores	Pain Score per Area		Drafting	10/05/18 6:33 PM
L51	Overall Relief Score			Drafting	08/28/18 10:53 AM
L52	No Data			Drafting	08/28/18 3:53 PM
L53	Pain Area			Drafting	08/28/18 3:57 PM
L54	Score			Drafting	08/28/18 3:57 PM
L55	Baseline			Drafting	08/28/18 3:57 PM
L56	HF10 Off			Drafting	08/28/18 4:37 PM
L57	HF10 On			Drafting	08/28/18 4:37 PM

Content Tag	Original Content	Latest Content	Comments	Status	Date Modified	
L58	Active Favorite			Drafting	09/03/18 5:40 PM	
L59	Add Favorite			Drafting	09/03/18 6:03 PM	
L60	Add Current Program			Drafting	09/03/18 6:03 PM	
L61	Add from Program History			Drafting	09/03/18 6:03 PM	
L62	Edit Favorite			Drafting	09/03/18 6:44 PM	
L63	Delete Favorite			Drafting	09/03/18 7:01 PM	
L64	Unable to Add Favorite			Drafting	09/03/18 7:01 PM	
L65	Setting			Drafting	10/08/18 1:40 PM	
L66	"[X]" Favorite		Where [X] = Favorite Name that's currently running for either Pathfinder or PTA is displayed in quotes, followed by the word Favorite. e.g., "Walking" Favorite; "Running" Favorite; "Sleeping" Favorite	Drafting	10/09/18 1:37 PM	
L67	HF10 is currently ON.			Drafting	10/09/18 12:46 PM	
L68	Favorite in Use			Drafting	10/10/18 10:58 PM	
L69	[X] Optimized		Where [X] = Pathfinder Program Name. e.g, "JD2 Optimized"	Drafting	10/10/18 10:58 PM	
L70	Connection Lost			Drafting	10/10/18 10:58 PM	
L71	Status		Label at top of home screen cards	Drafting	10/10/18 10:58 PM	
T1	EULA Agreement Copy		Currently with Nevro Legal	Drafting	07/17/18 12:51 PM	
T2	To use this app for controlling your HF10 Implant, you must agree to the End User License Agreement on the screen. If you are certain you wish to decline acceptance of these terms, please quit and uninstall the app from your phone.		Text for user who declines EULA - Nevro needs to review and determine if it is necessary	Drafting	07/17/18 12:51 PM	
T3	Your phone isn't compatible with HF10. You can recheck your phone or seek Support / Help.			Drafting	07/17/18 12:51 PM	
T4	Scanning for HF10 Implant...			Drafting	07/17/18 12:51 PM	
T5	Found HF10 Implant	HF10 Implant Found	Header	Drafting	07/31/18 11:16 AM	
T6	Tap on your HF10 Implant below to pair, ensuring the serial number matches the one on the card given to you after surgery.	Confirm the serial number matches your patient ID, then tap the PAIR button below	Instructional Text	Drafting	07/31/18 11:16 AM	
T7	Don't see your device? Tap below to start over.	Don't see your device?	Instructional Text	Drafting	07/31/18 11:16 AM	
T8	Attempting to Pair...		System Status	Drafting	07/17/18 12:51 PM	
T9	Congratulations!		Header	Drafting	07/17/18 12:51 PM	
T10	Your HF10 Implant has been paired successfully.		System Status	Drafting	07/17/18 12:51 PM	
T11	"Be sure your HF10 Implant is in "Discoverable" Mode. It will stay in this mode for two minutes. Tap below to start over and try again.		Instructional Text	Drafting	07/17/18 12:51 PM	
T12	We were unable to pair with your HF10 Implant because it is no longer in "Discoverable" mode. It will stay in this mode for two minutes. Tap below to start over and try again.	We were unable to pair with your HF10 Implant because it is no longer in "Discoverable" mode. It will stay in this mode for two minutes.	Instructional Text	Drafting	07/31/18 11:29 AM	
T13	Please provide your first name		Instructional Text	Drafting	07/17/18 12:51 PM	
T14	mm/dd/yyyy, hh:mm AM/PM		Date and Time for Last Data Sync (L11)	Drafting	10/08/18 2:05 PM	
T15	Your changes will be applied, but stimulation will remain OFF until you turn it ON.		Text for user making program level change while stim is off. GH: No longer relevant as of 10/5/2018.	Drafting	10/05/18 12:30 PM	
T16	Slide to Turn Off		Text instructions for turning stimulation off. No longer relevant in new UI.	Drafting	10/08/18 3:07 PM	
T17	Slide to Turn On		Text instructions for turning stimulation on. No longer relevant in new UI.	Drafting	10/08/18 3:07 PM	
T18	Your HF10 is setup to target the pain area(s) above	Your HF10 is set up to target the pain area(s) above	Instructional Text	Drafting	07/31/18 11:45 AM	
T19	Select activities you'd like to do.	Select activities you'd enjoy	Instructional Text	Drafting	07/31/18 11:52 AM	
T20	Allow push notifications so your HF10 can notify you with low battery alerts and survey reminders.		Value prop statement for allowing notifications.	Drafting	07/17/18 12:51 PM	
T21	To turn Notifications on later, go to Menu > Notification Settings and select the options you'd like. You can also enable them from your iOS Settings > Notifications and scrolling to HF10.		Help Topic Content	Drafting	07/17/18 12:51 PM	
T22	To prevent accidental changes to your therapy settings, you can choose to setup Touch ID or PIN		Instruction Text	Drafting	07/17/18 12:51 PM	
T23	Your PIN has been saved		Confirmation Text	Drafting	07/17/18 12:51 PM	
T24	Your phone is now setup for use with your HF10 Implant		Confirmation Text	Drafting	07/17/18 12:51 PM	
T25	Please recharge soon.			Drafting	10/05/18 12:53 PM	
T26	Your HF10 battery is at 10%. Please recharge soon.		GH: Removed 10/5 since 25, 10 and 5% warnings are being modified to be "Low" and "Depleted"	Drafting	10/08/18 3:11 PM	
T27	Please recharge now.			Drafting	10/05/18 1:34 PM	
T28	Your HF10 battery is now fully charged.			Drafting	07/19/18 3:46 PM	
T29	There is a problem connecting to your HF10.			Drafting	07/19/18 3:43 PM	
T30	What is your overall pain relief today?			Drafting	07/27/18 4:18 PM	
T31	What is your X pain score today?		The question is asked n number of times, with each time X equal to the respective pain area name from T32 array below.	Drafting	07/30/18 2:23 AM	

Content Tag	Original Content	Latest Content	Comments	Status	Date Modified	
T32	left foot, right foot, bilateral feet, left leg, right leg, bilateral leg, perineum, left buttock, right buttock, bilateral buttock, lower back, mid back, upper back, neck, left hand, right hand, left arm, right arm, bilateral arm, left shoulder, right shoulder, bilateral shoulder, groin, abdomen, chest, face		An array containing all 26 distinct pain area names for HF10 [0-25].	Drafting	07/30/18 2:23 AM	
T33	With HF10, your activity level has:			Drafting	07/30/18 11:12 AM	
T34	With HF10, your pain medication has:			Drafting	07/30/18 11:20 AM	
T35	Based on your answers, there is a recommended setting you should try.			Drafting	07/30/18 11:53 AM	
T36	Based on your answers, it's recommended you stay on our current settings. Your responses are correlated with data from your HF10 Implant, and you can view this in My Progress.			Drafting	07/30/18 11:53 AM	
T37	Based on your answers, it's recommended you test this program for another 48 hours. Your responses are correlated with data from your HF10 Implant, and you can view this in My Progress.			Drafting	07/30/18 11:53 AM	
T38	Your responses are correlated with data from your HF10 Implant, and you can view this in My Progress.			Drafting	07/30/18 12:06 PM	
T39	To maximize the benefit you get from HF10, we encourage you to check out more information in Tips & Guidance.			Drafting	07/30/18 12:06 PM	
T40	Are you sure you wish to exit the survey? Your responses won't be saved.			Drafting	07/30/18 12:34 PM	
T41	48 hours is recommended when trying out new program settings.	Staying with each program for 48 hours is recommended when trying out new program settings.		Drafting	07/31/18 12:17 PM	
T42	Take the short survey so we can provide recommendations to optimize your HF10.			Drafting	07/30/18 6:15 PM	
T43	Take the short survey so you can log your progress with HF10.			Drafting	07/30/18 6:15 PM	
T44	We encourage you to continue testing the current program for another X hours.			Drafting	08/28/18 11:34 AM	
T45	Turn on HF10 and begin taking surveys to see this data mapped above.			Drafting	08/28/18 3:53 PM	
T46	Turn on HF10 and begin taking surveys to see your results graphed over time.			Drafting	08/28/18 4:37 PM	
T47	No survey was taken, so we don't know what % pain relief you experienced during this program.			Drafting	08/28/18 4:37 PM	
T48	You can save up to 10 favorites			Drafting	09/03/18 5:12 PM	
T49	You can save [X] more favorites		Where X is the number of remaining favorite slots (10 total) that are available.	Drafting	09/03/18 5:14 PM	
T50	You can only save 10 favorites. Please delete one before adding another.			Drafting	09/03/18 5:14 PM	
T51	Choose a custom name (12 character max) and pick a color to label your favorite.			Drafting	09/03/18 6:08 PM	
T52	Are you sure you want to delete [X]?		Where X is the name of the Favorite the patient is wanting to delete.	Drafting	09/03/18 7:01 PM	
T53	1. Unplug the charger		Frame 1 in Pairing Animation	Drafting	10/07/18 6:41 PM	
T54	2. Press the power button until you hear a beeping sound		Frames 2 & 3 in Pairing Animation	Drafting	10/07/18 6:41 PM	
T55	3. Place charging coil over implant site on body		Frame 4 in Pairing Animation	Drafting	10/07/18 6:41 PM	
T56	4. When aligned, chirping will stop and battery indicator will light up		Frames 5 & 6 in Pairing Animation	Drafting	10/07/18 6:41 PM	
T57	Tap to turn HF10 ON			Drafting	10/09/18 12:41 PM	
T58	Your HF10 has been turned OFF.			Drafting	10/09/18 12:46 PM	
T59	Your HF10 has been turned ON.			Drafting	10/09/18 12:46 PM	
A1	Animation: How to Discover & Pair Your HF10 Implant		Looping animation depicting the steps for putting Implant into Discoverable Mode. See animation_discover-pair[1-10] in Zeplin for animation specs	Drafting	07/19/18 7:42 PM	
A2	Animation: HF10 Pairing Attempt		Animation Specs in Progress	Drafting	07/19/18 7:42 PM	
I1	Illustration: HF10 with checkmark		Placeholder image, use illustration in Zeplin comp as working placeholder	Drafting	07/19/18 8:04 PM	
I2	Illustration: HF10 with question mark icon		Placeholder image, use illustration in Zeplin comp as working placeholder	Drafting	07/19/18 8:04 PM	
I3	Illustration: HF10 with warning icon		Placeholder image, use illustration in Zeplin comp as working placeholder	Drafting	07/19/18 8:04 PM	
I4	Illustration: Generic type of activity depicted.		Placeholder image, use illustration in Zeplin comp as working placeholder	Drafting	07/19/18 8:04 PM	
I5	Illustration: iOS Phone with HF10 Notification Banner. Text is: It's Time to Optimize your HF10. Take the Survey.		Placeholder image, use illustration in Zeplin comp as working placeholder. Illustration for what a notification looks like with HF10.	Drafting	07/19/18 8:04 PM	

Content Tag	Original Content	Latest Content	Comments	Status	Date Modified	
I6	<i>Illustration: Why Allow Notifications</i>		Placeholder image, use illustration in Zeplin comp as working placeholder	Drafting	07/19/18 8:04 PM	
I7	<i>Illustration: Touch ID or other graphic representing authentication on iOS / Android.</i>		Placeholder image, use illustration in Zeplin comp as working placeholder	Drafting	07/19/18 8:04 PM	
I8	<i>Illustration: Trophy cup</i>		Placeholder image, use illustration in Zeplin comp as working placeholder	Drafting	07/19/18 8:04 PM	
I9	Illustration: Posterior view of gender-neutral silhouette of human body upon which pain area overlays will be placed.		See iOS/Android UI Kit in Zeplin	Drafting	07/19/18 7:52 PM	
I10	Illustration: Pain area overlay that will be used to visually indicate a pain area.		See iOS/Android UI Kit in Zeplin	Drafting	07/19/18 7:52 PM	
I11	Illustration I9 + Anterior view of same silhouette and placed next to it when pain areas from both front AND back are being targeted.		See iOS/Android UI Kit in Zeplin	Drafting	07/19/18 7:52 PM	
I12	Illustration: Activity Level Lifestyle Image for Survey Question		-	Drafting	07/30/18 11:58 AM	
I13	Illustration: Pain Medication Lifestyle Image for Survey Question		-	Drafting	07/30/18 11:58 AM	
I14	Illustration: Pathfinder Result 1		-	Drafting	07/30/18 12:28 PM	
I15	Illustration: Pathfinder Result 2		-	Drafting	07/30/18 12:28 PM	
I16	Swipe Left-Right Gesture		See iOS/Android UI Kit in Zeplin	Drafting	07/19/18 7:55 PM	
I17	Program Adjustment		See iOS/Android UI Kit in Zeplin. Contains #s 1-5, representing each of the five programs programmed for the patient.	Drafting	07/19/18 7:55 PM	
I18	Illustration: Pathfinder Result 3			Drafting	07/30/18 12:28 PM	
I19	Illustration: PTA Result 3			Drafting	07/30/18 12:28 PM	
I20	Illustration: PTA Result 3			Drafting	07/30/18 12:34 PM	
I21	Illustration: Pathfinder Result 3			Drafting	07/31/18 10:41 AM	
F1	Form Field		Empty by default	Drafting	07/19/18 7:55 PM	
F2	Form Field		Appears after tapping B20	Drafting	07/19/18 7:55 PM	
N1	Home		App Navigation	Drafting	07/18/18 3:00 PM	
N2	My Progress		App Navigation	Drafting	07/18/18 3:00 PM	
N3	Favorites		App Navigation	Drafting	07/18/18 3:00 PM	
N4	Menu		App Navigation	Drafting	07/18/18 3:00 PM	